

PRODUCT SPECIFICATION

Spec code: NOF-C0940
Product name: Semperguard® Style
Date of issue: 20/05/14
Article number: 8167 6413.
 . = 3 (size S), 5 (size M), 7 (size L), 9 (size XL)

GENERAL INFORMATION

Type: single-use disposable examination or protective glove
Marking: all information on dispenser box
Shape: ambidextrous – straight fingers
Material: Nitrile Butadiene Rubber (NBR)
Colour: black
Inner: online chlorinated/powder free
Outer: no treatment
Cuff/Surface: rolled cuff/fingertip textured

Thickness: (double measured)	Finger: min. 0.20 mm	Palm: min. 0.14 mm	Cuff: min. 0.10 mm
--	-----------------------------	---------------------------	---------------------------

Residual powder: [according to EN 455-3]	≤ 2 mg / glove
--	----------------

Force at break:	med. 6.0 N (during shelf life)
------------------------	-----------------------------------

Length:	median 240 mm
----------------	---------------

Available Sizes:	S (6-7)	M (7-8)	L (8-9)	XL (9-10)
-------------------------	---------	---------	---------	-----------

Storage: 3 years storage life (date of expiration see on the box). Keep area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper-ions discolour the glove. Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.

STANDARDS, GUIDELINES AND QUALITY CERTIFICATION

Quality certification: ISO 9001, ISO 13485
Conformity to guidelines: Personal Protective Equipment 89/686/EEC: CE cat. III
 Medical Device Directive 93/42/EEC: CE Class I

Conformity to standards:	EN 420 EN 455 ASTM D6319	EN 374-2: specified level 2 - AQL 1.5 level 3: AQL <0.65 (G1) level 2: AQL <1.5 (G1) level 1: AQL <4.0 (S4)	EN 374-3: Class 6 Sodium Hydroxide 40%
---------------------------------	--------------------------------	---	---

Conformity to food regulation: 1935/2004






WARNING AND REPORTING SYSTEM

Potentially allergenic ingredients: Dithiocarbamate type, Zinc-mercaptobenzothiazol
 (For further information, the actual version of the ingredients list is given on request.)

Medical Device Vigilance and Reporting System: According to the official reporting criteria of MDD, incidents caused by examination gloves or surgical gloves must be reported immediately to our medical device reporting officer. Tel.: +43 2630 310 ext. -510; Fax: -549

Caution! This product contains accelerators which may cause allergic reactions!

Established (PM/DM):	Reviewed (QA/DM):	Approved (HoM/DM):	Customer approval
 L. Millonig, MA	 Dr. R. Stempfer	 Mag. B. Schnetzlinger	