

# PRODUCT SPECIFICATION

**Spec code:** NOT-C0367 290mm  
**Product name:** Semperguard® Nitrile Xtension  
**Date of issue:** 13/10/14  
**Article number:** 8167 8086.  
 . = 3 (size S), 5 (size M), 7 (size L), 9 (size XL)

## GENERAL INFORMATION

**Type:** single-use disposable examination or protective glove  
**Marking:** all information on dispenser box  
**Shape:** ambidextrous – straight fingers  
**Material:** Nitrile Butadiene Rubber (NBR)  
**Colour:** blue  
**Inner:** online chlorinated/powder free  
**Outer:** no treatment  
**Cuff/Surface:** rolled cuff/fully textured

<b>Thickness:</b> <small>(double measured)</small>	<b>Finger:</b> min. 0.18 mm	<b>Palm:</b> min. 0.16 mm	<b>Cuff:</b> min. 0.10 mm
---	-----------------------------	---------------------------	---------------------------

<b>Residual powder:</b>	≤ 2 mg / glove
-------------------------	----------------

<b>Force at break:</b> <small>[according to EN 455-2]</small>	med. 6.0 N <small>(during shelf life)</small>
--	--

<b>Length:</b>	median 290 mm
----------------	---------------

<b>Available Sizes:</b>	S (6-7)	M (7-8)	L (8-9)	XL (9-10)
-------------------------	---------	---------	---------	-----------

**Storage:** 3 years storage life (date of expiration see on the box). Keep area cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper-ions discolour the glove. Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.

## STANDARDS, GUIDELINES AND QUALITY CERTIFICATION

**Quality certification:** ISO 9001, ISO 13485  
**Conformity to guidelines:** Personal Protective Equipment 89/686/EEC: CE cat. III  
 Medical Device Directive 93/42/EEC: CE Class I

<b>Conformity to standards:</b>	EN 420	EN 374-2: specified	EN 374-3: Class 6
	EN 455	level 2 - AQL 1.5	Sodium Hydroxide 40%
	ASTM D6319	level 3: AQL <0.65 (G1) level 2: AQL <1.5 (G1) level 1: AQL <4.0 (S4)	




**Conformity to food regulation:** 1935/2004 

## WARNING AND REPORTING SYSTEM

**Potentially allergenic ingredients:** Dithiocarbamate type  
*(For further information, the actual version of the ingredients list is given on request.)*

**Medical Device Vigilance and Reporting System:** According to the official reporting criteria of MDD, incidents caused by examination gloves or surgical gloves must be reported immediately to our medical device reporting officer. Tel.: +43 2630 310 ext. -510; Fax: -549

**Caution!** This product contains accelerators which may cause allergic reactions!

<b>Established (PM):</b>  L. Millonig, MA	<b>Reviewed (QA):</b>  Dr. R. Schaller	<b>Approved (R&amp;D):</b>  Dr. A. Weinert
--	---	--